



1250 Sheppard Ave. West
Toronto, ON M3K 2A6

T 1.800.387.3713
F 416.661.3447

order@waynesafety.com

www.waynesafety.com

SPECIFICATION FOR 15I-750PF NITRILE EXAMINATION GLOVE

1.0 Product:

Purple Nitrile, Powder Free, Textured, 3mil, Medical Grade For Examination Use.

2.0 Physical Dimensions:

Size	Length (mm)	Width (mm)	Thickness (mm) (min)	
			Palm	Finger
S	230mm±10mm	80±10	0.06	0.07
M	230mm±10mm	95±10	0.06	0.07
L	230mm±10mm	110±10	0.06	0.07
XL	230mm±10mm	120±10	0.06	0.07
XS	300mm±10mm	70±10	0.06	0.07

3.0 Physical Requirements:

Before Aging:

Tensile Strength: 14 MPA; Ultimate Elongation: 500% (min)

After Accelerated Aging:

Tensile Strength: 14 MPA; Ultimate Elongation: 400% (min)

Note: For gloves that are older than 6 months from the date of manufacture or for which the date of manufacture cannot be determined, no accelerated aging shall be performed. If such gloves are tested, the tensile strength and ultimate elongation shall not be less than the "After Accelerated Aging" values.

4.0 Freedom from holes:

According to ISO 2859, G-II, multiple sampling plan for normal inspection, AQL 1.5

5.0 Powder Residues: Not more than 2mg per glove.

Conforms to EN375 and EN455-2 Standards
